
# Trooper Count Scale Instruction Manual

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**Declaration of Conformity** We, Ohaus Corporation, declare under our sole responsibility that the balance models listed below marked with "CE" - are in conformity with the directives and standards mentioned.

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Date: March 28, 2003

Ted Xia  
President  
Ohaus Corporation  
Pine Brook, NJ USA

Johan Dierbach  
General Manager  
Ohaus Europe  
Greifensee, Switzerland

**ISO 9001 Registration for Ohaus Corporation.** Ohaus Corporation, USA, was examined and evaluated in 1994 by the Bureau Veritas Quality International, BVQI, and was awarded ISO 9001 registration. This certifies that Ohaus Corporation, USA, has a quality system that conforms with the international standards for quality management and quality assurance (ISO 9000 series). Repeat audits are carried out by BVQI at intervals to check that the quality system is operated in the proper manner.

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**NOTE:** THIS EQUIPMENT HAS BEEN TESTED AND FOUND TO COMPLY WITH THE LIMITS FOR A CLASS A DIGITAL DEVICE, PURSUANT TO PART 15 OF THE FCC RULES.

THESE LIMITS ARE DESIGNED TO PROVIDE REASONABLE PROTECTION AGAINST HARMFUL INTERFERENCE WHEN THE EQUIPMENT IS OPERATED IN A COMMERCIAL ENVIRONMENT. THIS EQUIPMENT GENERATES, USES, AND CAN RADIATE RADIO FREQUENCY ENERGY AND, IF NOT INSTALLED AND USED IN ACCORDANCE WITH THE INSTRUCTION MANUAL, MAY CAUSE HARMFUL INTERFERENCE TO RADIO COMMUNICATIONS. OPERATION OF THIS EQUIPMENT IN A RESIDENTIAL AREA IS LIKELY TO CAUSE HARMFUL INTERFERENCE IN WHICH CASE THE USER WILL BE REQUIRED TO CORRECT THE INTERFERENCE AT HIS OWN EXPENSE.

THIS DIGITAL APPARATUS DOES NOT EXCEED THE CLASS A LIMITS FOR RADIO NOISE EMISSIONS FROM DIGITAL APPARATUS AS SET OUT IN THE INTERFERENCE-CAUSING EQUIPMENT STANDARD ENTITLED "DIGITAL APPARATUS", ICES-003 OF THE DEPARTMENT OF COMMUNICATIONS.

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**Unauthorized changes or modifications to this equipment are not permitted.**

# TABLE OF CONTENTS

1. GETTING TO KNOW YOUR TROOPER COUNT SCALE .....	3
1.1 Introduction .....	3
1.1.1 Features .....	3
2. INSTALLATION .....	3
2.1 Unpacking and Checking .....	3
2.2 Selecting the Location .....	4
2.3 Connecting Power .....	4
2.3.1 Battery Installation .....	4
2.3.2 AC Power .....	4
2.4 Leveling the Scale .....	4
2.5 Stabilization .....	4
3. OPERATION .....	5
3.1 Overview of Controls and Display Functions .....	5
3.2 Basic Weighing .....	7
3.2.1 Turning On the Scale .....	7
3.2.2 Turning Off the Scale .....	7
3.2.3 Zero Operation .....	7
3.2.4 Weighing .....	7
3.2.5 Manual Tare Operation .....	8
3.2.6 Pre-Set Tare Operation .....	8
3.2.7 Gross/Tare Recall Operation .....	8
3.2.8 Unit Switch Operation .....	8
3.3 Basic Counting .....	9
3.3.1 Positive Sampling .....	9
3.3.2 Negative Sampling .....	9
3.3.3 Positive Counting .....	10
3.3.4 Negative Counting .....	10
3.3.5 Using Keypad to Edit Counting Data .....	11
3.4 Library Mode .....	11
3.4.1 Library Operation .....	11
3.4.2 Entering Data Directly into Library .....	11
3.4.3 Saving Library Data .....	13
3.4.4 Saving Library Data in Another Location .....	13
3.4.5 Using the Library .....	14
3.5 Auto Optimization .....	14
3.6 Auto Add Mode .....	14

## TABLE OF CONTENTS (Cont.)

3.7 RS232 Communications .....	15
3.7.1 RS232 User Command Table .....	15
3.7.2 Output Formats .....	15
3.7.3 RS-232 Pin Out .....	15
4. SETUP.....	16
4.1 Setup Protection .....	16
4.2 Menu Structure .....	16
4.3 Control Functions .....	17
4.3.1 Menu Operation .....	17
4.4 Checking the Menu Settings .....	17
4.5 Setup Menu .....	18
4.6 Read Menu .....	19
4.7 Print Menu .....	19
4.8 Count Menu .....	21
4.9 Lockout Menu .....	22
5. CALIBRATION AND SEALING .....	23
5.1 Span Calibration .....	24
5.2 Linearity Calibration .....	25
5.3 Legal for Trade Operation and Sealing .....	27
6. CARE AND MAINTENANCE .....	28
6.1 Troubleshooting .....	28
6.2 Error Codes List .....	29
6.3 Service Information .....	29
6.4 Accessories .....	29
7. TECHNICAL DATA .....	28
7.1 Counting Limits .....	28
7.2 Specifications .....	30

## **1. GETTING TO KNOW YOUR TROOPER COUNT SCALE**

### **1.1 Introduction**

Thank you for deciding to purchase a Trooper Count Scale from Ohaus. The Ohaus Trooper Count Scale is a rugged, reliable, electronic counting/weighing scale designed for easy operation.

Behind your instrument stands Ohaus, a leading manufacturer of precision Indicators, Scales and Balances. An After-market Department with trained instrument technicians is dedicated to providing you with the fastest service possible in the event your instrument requires servicing. Ohaus also has a Customer Service Department to answer any inquiries regarding applications and accessories.

To ensure you make full use of the possibilities offered by your Trooper Count Scale, please read the manual completely before installation and operation.

#### **1.1.1 Features**

Major features include:

- Three backlit LCDs
- Durable cast aluminum internal frame, stainless steel weighing pan, ABS housing
- 21 button operation with numeric keypad
- Flexible unit switching-lb/kg/oz/g
- RS232 serial communication interface
- Battery operation with internal 6 C-cell batteries
- AC power adapter
- Auto shut off for power saving
- Library feature saves counting data for 60 items. (Part number, APW, tare values, over/under values, and count totals.)
- Positive or negative sampling and counting

## **2. INSTALLATION**

### **2.1 Unpacking and Checking**

Open the package and remove the instrument and the accessories. Check the instrument for transport damage. Immediately inform your Ohaus dealer if you have complaints or if parts are missing. Your Trooper Count Scale package should contain:

- Trooper Count Scale
- Weighing Platform
- AC power adapter
- Instruction Manual
- Warranty card
- Lead seal kit and capacity labels

Store all parts of the packaging. This packaging ensures the best possible protection for the transport of your instrument.




**NOTE:** Remove the two shipping screws if present on top of the unit.

# Trooper Count Scale

---

## 2.2 Selecting the Location

The Scale should be used in an environment which is free from dust, water, corrosives, vibration or temperature extremes as these factors can affect scale accuracy. The Scale should be located on a stable level surface and kept away from vibrating sources such as large machinery.

## 2.3 Connecting Power

### 2.3.1 Battery Installation (If desired)

Remove battery cover and insert 6 Alkaline C-type batteries into the battery holder making sure that the batteries are properly orientated. The negative battery terminals should face the spring contacts.

**NOTE:** It is recommended that when the Trooper Count Scale is operated from batteries, the Auto-Off Timer feature be turned on to extend battery life. When setting up the Scale, refer to the Setup menu, section 4.5.


### 2.3.2 AC Power

Verify that the AC power source matches the AC adapter rating. Connect the Adapter power cord at the rear of the Scale and plug the Adapter into a convenient power outlet.

## 2.4 Leveling the Scale

To compensate for small irregularities or inclinations at the location, the Scale can be leveled. The Scale is equipped with a level indicator located at the rear. Adjust the leveling feet until the air bubble in the indicator is centered as shown.

**NOTE:** The instrument should be leveled each time its location is changed.


## 2.5 Stabilization

Before initially using the Scale, allow time for it to adjust to its new environment. The recommended warm up period is five (5) minutes if the Scale is already stabilized to room temperature.

### 3. OPERATION

#### 3.1 Overview of Controls and Display Functions

No.	Designation	Function
1	Weight Display	<p>Displays current weight on scale.</p> <p>Indicators identify active weight unit and Net or Gross weight:  <b>kg</b> (Kilogram) – Unit of measure.  <b>g</b> (Gram) - Unit of measure.  <b>lb</b> (Avoirdupois Pound) - Unit of measure.  <b>oz</b> (Avoirdupois Ounce) - Unit of measure.  <b>Net</b> (Net weight) – Gross weight minus Tare.</p>
2	Data Display	<p>Displays user selected data.</p> <p>Indicators show active data mode:  <b>B/G</b> (Brutto/Gross) – Total weight of items.  <b>%AC</b> (Percent Count accuracy) – Estimated accuracy of the count based on sample size and scale resolution.  <b>T</b> (Tare value) – Weight of the container .  <b>PT</b> (Preset Tare value) – Tare value recalled from memory or entered via keypad.  <b>APW</b> (Average Piece Weight) – Average weight of one item.  <b>MR</b> (Memory Recall) – Accumulated quantity stored in memory.</p>
3	Count Display	<p>Displays current count of items on the scale.</p> <p>Indicators show count status:  <b>M</b> (Memory) – Accumulated quantity in memory.  <b>&lt;</b> (Under) – Current quantity is less than lower limit.  <b>&gt;</b> (Over) – Current quantity is more than upper limit.  <b>&lt; &gt;</b> (Target range) – Current quantity is between upper and lower limits.  <b>PCS</b> (Pieces) – A count value is being shown.</p>

## 3.1 Overview of Controls and Display Functions (Cont.)

No.	Designation	Function
4	Numeric Keypad	Enter specific values for tare, APW and other data entries.
5	Data button	Selects between data modes in the Data display. In edit mode, serves as a "Go Back" or backspace key.
6	Clear button	During editing, clears the blinking entry.
7	M+ button	Adds the displayed quantity into the accumulating memory. During editing, selects the next higher option or number.
8	M- button	Subtracts the keyboard entry from the accumulating memory. In sampling mode, enters negative sampling. During editing, selects the next lower option or number.
9	Enter button	Enters edit function of current mode. If current mode does not support edit functions, an error beep will occur. In edit mode, accepts entry.
10	TARE button	Tares the weight on the pan or accepts the keypad tare entries.
11	LIBRARY <i>Save</i> button	Short press displays Library ID (blinking) and the part number. Repeat press returns to weigh/count mode. Long press (3 seconds) saves library changes.
12	SAMPLE <i>Menu</i> button	Short press enters the sampling mode. Long press (3 seconds) enters scale setup menu. In menu mode, a short press will return to weigh mode.
13	PRINT <i>Units</i> button	Short press sends data to RS232 port. Long press (3 seconds) steps through active weigh units. Release on desired unit.
14	ON/ZERO <i>Off</i> button	Short press turns the scale on or zeros the display. Long press (3 seconds) turns the scale off.


## 3.2 Basic Weighing

The Trooper Count Scale has been factory calibrated and is ready for operation. Ohaus recommends that the Scale be re-calibrated before use to ensure optimal accuracy. You can operate the scale at this point using the factory default settings. You can check the default menu settings on page 14. All bolded items on the menu are the factory default settings. The scale will operate in the default mode. If you want to change the settings, continue with Section 4, Setup.

When the Scale is positioned for operation, follow the operational procedure outlined next.


### 3.2.1 Turning On the Scale

Press and hold **ON/ZERO Off** button until the LCD display appears. The display cycles as shown below displaying "Ohaus", with model number, serial number, software revision number and weighing mode display.


### 3.2.2 Turning Off the Scale

To turn the Scale off, press and hold the **ON/ZERO Off** button until OFF is displayed.

### 3.2.3 Zero Operation


Press the **ON/ZERO Off** button to zero the weight display. The gross weight is also shown in the DATA display.

**NOTE:** A pre-set tare value will not be zeroed.

### 3.2.4 Weighing


Place item to be weighed on the Scale. The display indicates a sample of 1.5kg, gross weight.


## 3.2.5 Manual Tare Operation

When weighing a sample that must be held in a container, taring stores the container weight memory.


Place the container on the Scale. Display shown is 0.5kg.

Press **TARE** button. The container is tared, the Net weight is shown in the WEIGHT display, and the Tare weight is shown in the DATA display.


To clear the tare with an empty pan, press the **TARE** button.


## 3.2.6 Pre-Set Tare Operation

Pre-Set Tare is a tare value that is entered using the numeric keypad. To use Pre-Set Tare, repeatedly press the **Data** button until the PT indicator is shown in the DATA display. The DATA display indicates OFF.


Press the **Enter** button to toggle the Pre-Set Tare on or off. The current Pre-Set Tare value (if any) is displayed. Enter a numeric value using the keypad, (the COUNT display is blinking), then press the **Enter** button. The Net weight is shown in the WEIGHT display and the tare weight is shown in the DATA display. Display shown is 0.5kg.

### Alternate Method:

Pre-Set Tare can also be entered using the numeric keypad, then press the **TARE** button.

**NOTE:** The Pre-Set Tare value is overwritten when a new value is entered.


## 3.2.7 Gross/Tare Recall Operation

Repeatedly press the **Data** button to cycle through Gross and Tare readings.

## 3.2.8 Unit Switch Operation

Press and hold **PRINT Units** button until display changes to selected measuring unit. Display indicators at the bottom of the WEIGHT display indicate the active weighing unit. Depending on which units are enabled in the Read menu, you have a choice of g, lb, kg or oz. The display sample indicates a 1kg load changed to lb.

### 3.3 Basic Counting

Read this procedure entirely before counting pieces. Piece sampling and counting can be further enhanced using menu selections in the Count menu. There are several methods that are available for sampling.



- Positive Sampling - given number of pieces are placed on the scale, the weight is divided by the number of pieces to calculate the Average Piece Weight (APW). Additional pieces are then placed on the scale and counted and displayed.
- Negative Sampling - calculates the APW when a given number of pieces are removed from the total quantity.
- The APW can also be entered through the numeric keypad. As parts are placed on the scale, the count is immediately displayed.
- The Trooper Count will support positive counting or negative counting.

#### 3.3.1 Positive Sampling

If a weighing container is to be used:

- If the tare value is known, enter a Pre-Set Tare.
- If the tare value is not known, place the container on the pan and press the **TARE** button.

Press the **SAMPLE Menu** button.

The display indicates to place 10 samples on the scale. Place the desired number of sample pieces in the container. If you have placed 10 pieces in the container, press the **Enter** button, **or** by pressing the numeric keys you may directly enter an alternate sample size. Increased sample size results in better results. If you make a mistake, press the **Clear** button and start again. Press the **Enter** button. The example indicates a sample size of 10 pieces.

#### NOTES:

1. If **APW** is less than minimum, [Err3] [LOW] [APW] is displayed. See section 6.2.
2. If **Ref. Wt.** is less than minimum, [Err4] [LOW] [REF] is displayed. See section 6.2.

#### 3.3.2 Negative Sampling


Negative sampling permits the determination of the total quantity on the pan by removing a known number of sample pieces. The total quantity may be placed directly on the Scale or in a container. The Trooper Count Scale automatically calculates the remaining pieces after the sample amount is removed. When a container is used, you may enter a preset tare weight or simply tare the empty container before starting.

## 3.3.2 Negative Sampling (Cont.)


If a weighing container is to be used:

- If the tare value is known, enter a Pre-Set Tare.
- If the tare value is not known, place the container on the pan and press the **TARE** button.

Press the **SAMPLE Menu** button.

Press the **M-** button.


Place the total number of pieces on the pan, then press the **Enter** button.

Remove 10 pieces as instructed, or use the keypad to change the number of pieces to be removed, then press **Enter** button.

Place the pieces that were removed back on the scale to display the total number of pieces.

**NOTE:** If AUTO-ADD is ON, the Scale may direct you to remove additional pieces if the reference weight is too small.

## 3.3.3 Positive Counting

After sampling, place additional pieces on the scale and view the piece count in the COUNT display.

## 3.3.4 Negative Counting

After sampling, place the entire quantity of pieces on the scale. Press the **ON/ZERO Off** button. As pieces are removed, the quantity removed is shown in the COUNT display as a negative number.

### 3.3.5 Using Keypad to Edit Counting Data

The APW, Pre-Set Tare (PT) and the accumulated count (MR) can be viewed and edited in the DATA window. While in weigh mode repeatedly press the **Data** button until the APW indicator is lit. A new APW value may be entered via the keypad and accepted with an **Enter** press. PT and MR can also be similarly edited.

### 3.4 Library Mode

The Trooper Count scale allows the storage of part/counting data in 60 locations, designated ID 00 to ID 59. Each location includes: part number, Average Piece Weight (APW), a pre-set tare value, over and under count values and an accumulated piece count.

Data in the library is saved in a manner similar to a computer file. Temporary editing of the data will not affect the saved data. When desired, the temporary edits can be saved, allowing the scale to be turned off without losing the changes. Data in ID 00 may not be saved permanently but the data can be easily copied to another location.

The Library, ID 01 to ID 59, can be locked (see 4.9 Lockout Menu) to prevent unauthorized changes to the data. The accumulated count will not be locked allowing it to be updated. When turned on, the scale defaults to ID 00.

#### 3.4.1 Library Operation

To enter the Library, momentarily press the **LIBRARY Save** button. A short press in Library mode will exit to weigh mode.


- **M+** and **M-** buttons advance through the blinking choices.
- The **Enter** button accepts the blinking choice.
- The **Data** button returns to the previous choice.
- The **Clear** button will clear blinking data or back space during data entry.
- Three rapid beeps indicates an error, an improper button press or data value.

#### 3.4.2 Entering Data Directly into Library


The following example is used to illustrate how to enter data in a particular library location.

Short press on the **LIBRARY Save** button.

Short presses on the **M+** or **M-** buttons will display library locations which are already used.


Using the numeric keys, enter the desired library number. The sample shown is an empty location, ID 22. Press the **Enter** button to accept ID 22.


Press the **Enter** button, to edit the part number or press **M+** to skip.

## 3.4.2 Entering Data Directly into Library (Cont.)


Using the numeric keys 0-9, enter the desired part number. For this example, we will use 123000000. You will notice that the part number can be up to 12 digits and utilizes the DATA and COUNT display windows. Press the **Enter** button to confirm the entry.


Press the **M+** button, the display advances to APW. Press the **Enter** button to access APW value.


Using the numeric keypad, enter the APW. Sample shown indicates an average piece weight of 2g. Press the **Enter** button to accept.


Press the **M+** button to advance to setting the preset tare. Press the **Enter** button to access the tare value.


Using the keypad, enter the desired tare weight. The sample shown is 50g. Press the **Enter** button to accept.


Press the **M+** button to advance to Under count. Press **Enter** button to access the under count value.


Using the numeric keypad, enter the desired under count limit. For this sample, we will use 95. Press the **Enter** button to accept.

Press **M+** to advance to over count. In this sample, we are using 105 pieces.


Press **Enter** button to access over count value. Using the numerical keypad, enter the desired over count limit (must be greater than under count value). Press the **Enter** button to accept.


### 3.4.2 Entering Data Directly into Library (Cont.)

Press **M+** to advance to count. Count is the accumulated piece count in memory. Press **Enter** button to edit count value.

Using the numeric keypad, enter the desired count. The sample illustrates 250 pieces. Press the **Enter** button to accept.


To exit library editing, momentarily press the **Library**


### 3.4.3 Saving Library Data

Temporary changes to a library location need to be permanently saved or the changes will be lost when the scale is turned off. The scale will recommend saving temporary changes when the library location is changed or when the scale is turned off.


To save the settings, use a long press on the **Library Save** button. The library location (ID 22) and UPDATE YES is displayed.

Press **Enter** button to accept save location. Press **Enter** button again to save changes.


The UPDATE DONE display confirms a successful save.


### 3.4.4 Saving Library Data in Another Location

You have the option of copying all data from one library location into another library location. In the following example, we will save data from library location 22 to 52. Library 22 is shown.

Make a long press on **LIBRARY Save** button until display indicates UPDATE.


Enter the desired new location directory using the numeric keypad or by pressing the **M+** button to advance to the next open location. The sample is location 52. Press the **Enter** button.


## 3.4.4 Saving Library Data in Another Location (Cont.)


Press the **Enter** button to save data.

UPDATE DONE, then the ID number with IN USE will be displayed momentarily followed by the counting display.

At this point, all data from Library 22 has been copied into Library location 52.

## 3.4.5 Using the Library

To enter the Library, momentarily press the **LIBRARY Save** button. Use the keypad or the **M+** button to enter a location ID that you wish to use. Press the **LIBRARY Save** button again to select the new ID. ID XX IN USE will display indicating that the data stored in that memory is now active. If an APW and a Preset Tare was set, placing a sample filled container on the scale will give an immediate count.

If the Library is not locked, it is possible to edit the data in the current ID. The changes need not be saved if they were for temporary use. If the Library is locked, you must change to ID 00 for simple counting use.

## 3.5 Auto Optimization


Automatic reference optimization results in a more accurate count by increasing the reference weight without the need to count additional parts. A higher reference weight is important when there is a risk of inconsistent piece weights or if the reference weight is close to the minimum (See section 7.1).

Auto Optimization uses an initial APW to count additional pieces that are placed on the scale. After a few seconds, the message A-OPt appears in the Data display as the new higher reference weight is used to re-calculate the APW. The process can be repeated as long as the additional weight is less than the previous reference weight. Once this limit is exceeded Auto-Opt is turned off. This function can be activated or deactivated in the Count menu. The factory default setting is on.

## 3.6 Auto Add Mode

Auto Add Mode is setup for a more accurate count procedure by increasing the minimum reference weight and the minimum APW limits (see section 7.1). The scale requests additional pieces if the reference weight is less than the recommended Auto-Add minimum weight in Table 7-2. The instructions to add these additional pieces may be ignored. This function can be activated or deactivated in the Count menu. The factory default setting is off.

Place the stipulated number of pieces on the scale (e.g. 2 pieces). The scale then determines the new APW based on the total pieces

### 3.7 RS232 Communication

All communication is accomplished using standard ASCII format. Characters shown in the following table are acknowledged by the Scale. Invalid command response "ES" error indicates the Scale has not recognized the command. Commands sent to the Scale must be terminated with a carriage return (CR) or carriage return-line feed (CRLF). Data output by the Scale is always terminated with a carriage return-line feed (CRLF). See section 4.7 for RS232 and Print Setup.

#### 3.7.1 RS232 User Command Table

Command Character	Description (LFT is OFF)
<b>P</b>	Same as pressing PRINT button. Sends print data to peripheral device.
<b>V</b>	Output software version. For example: "TC30 Sr. 1.20"
<b>T</b>	Same as pressing TARE. If the Scale accepts the command, it will respond "OK".  If the Scale is in motion, it will output "CANT TARE" until the Scale is stable. When stable, tare is accepted. (NOTE: If LFT is ON, special tare rules apply.)
<b>Z</b>	Same as pressing ZERO button. If the Scale accepts the command, it will respond "OK".

**NOTE:** Any other command will produce the answer "ES".

#### 3.7.2 Output Formats

Data output can be initiated by: 1) By pressing **PRINT Units** button or, 2) Sending a print command ("P") from a computer. The two available output styles are shown below.

<pre>OHAUS TROOPER  COUNT TC3RS    SN: 222222  PART# = 123456789000 COUNT = 38 TOTAL = 188 APW = 5.000 g NET = 190 g TARE = 10 g GROSS = 200 g</pre>	<pre>OHAUS TROOPER  COUNT TC3RS    SERIAL No. : 222222  (WEIGHT UNIT = g) PART No.  COUNT    TOTAL  APW    NET    TARE    GROSS 123456789000  138    326    5.000  690    10     700</pre>
"Column" Print Style	"Line" Print Style

#### 3.7.3 RS-232 Pin Out.

The adjacent table illustrates the pin-out connections on the Female 9-pin D-sub connector.

1	N/C
2	Data In (RXD)
3	Data Out (TXD)
4	N/C
5	Ground
6, 7, 8, 9	N/C

## 4. SETUP


For first time setup, step through all menus and set the parameters as desired.

### 4.1 Setup Protection

The Trooper Count Scale is equipped with menus which permit certain functions to be locked out (not changed) during operation. If you intend to physically lock out changes to the setup selections, you will need to access the hardware Lock Switch located on the main circuit board following the setup procedure.

### 4.2 Menu Structure

The following table illustrates the menu structure in the Trooper Count Scale.

**NOTE:** Bolded items indicate factory defaults.

\* LFT on setting

### 4.3 Control Functions

The Scale has six menus; CAL, SETUP, READ, PRINT, COUNT and LOCSW which are entered by pressing and holding the **SAMPLE Menu** button until CAL MENU is displayed, then releasing it. During setup, the **M-**, **M+**, **Enter** and **Data** buttons are used to traverse and enter the menus.


#### 4.3.1 Menu Operation

- To enter menus, press and hold **SAMPLE Menu** button. A short press in menu mode will exit to the weigh mode.
- **M+** and **M-** buttons advances through blinking choices.
- **Enter** button accepts blinking choice.
- **Data** button returns to the next higher menu level. At top level, it exits the menu operation.
- Three rapid beeps indicates an improper button press.


### 4.4 Checking the Menu Settings

Each of the menus in the Trooper Count Scale can be quickly checked to see if any of the parameters need to be changed to match your requirements. The default settings are shown as bolded entries on the menu structure, section 4.2. Review the menu structure to determine what settings have to be changed.


To enter the menus, press and hold the **SAMPLE Menu** button until CAL MENU is displayed.


Repeated pressing of the **M-** button will traverse through the (6) menus in the order shown. You may go back by using the **M+** button. To leave the menus, press the **Data** button.


To enter a selected menu (blinking), press the **Enter** button, then scroll through all settings by repeatedly pressing the **M+** or **M-** buttons. In this manner, you can quickly check all settings for any particular menu.


All menu settings are saved once they are changed.

LOC ON indicates that the menu has been locked in the Loc Menu. You may view the menu parameters, but they may not be changed. See section 4.9 Lock-out menu for details.

The LOC MENU SW ON appears if the lockswitch underneath the scale is set on. This prevents the settings from being changed (as in a Legal For Trade application). See section 5.3 for details.

## 4.5 Setup Menu

### Procedure

With the Scale ON, press and hold the **SAMPLE Menu** button until CAL MENU is displayed.

**NOTE:** The calibration menu is covered in Section 5.

Press **M-** button to advance to SETUP MENU, then press **Enter** button. Navigate through the menu as explained in section 4.3.1.

Press **Enter** button to get into the RESET sub-menu. YES will change all menu settings to their default settings. Library data will not be affected.

ON will set the Scale to LFT configuration. In the ON setting, the zero setting will be limited to +/- 2%. In the OFF setting, the zero setting is 100%. Also, in the ON setting, the LFT readability is set. In the OFF position, the readability is set to the Enhanced value.










Selected calibration units are either lb or 1000g (kg).

A selection of 1/3, 2/3 or full span calibration values are available. The numeric values of the calibration points are displayed for the specific Scale capacity and calibration unit.

Three setting are available, ALL, SIGNAL and OFF. ALL = the beeper sounds when buttons are pressed or when error or target signals are needed. SIGNAL = beeps on errors (three beeps) or at count target (continuous). OFF = is always off.

Backlight has three available settings; ON, OFF, AUTO. ON = backlight is always on. AUTO = backlight turns on with any key press or a weight change (5d) and remains on for 5 seconds. OFF = backlight is always off. OFF or AUTO is recommended for battery operation.

ON = Scale turns off in 5 minutes if no key press or weight changes (5d) are detected. OFF = auto off feature is disabled. ON is recommended for battery operation.




## 4.6 Read Menu

### Procedure


If you have just finished changing settings in the SETUP menu, Press the **Data** button to go back to the SETUP MENU, SETUP should be blinking.

Press **M-** button to advance to the READ menu. In this menu, you are able to turn the measuring units on or off and set the Auto Zero tracking. Press the **Enter** button. Navigate through the menu as explained in section 4.3.1.

Grams (Unit g) - Select ON or OFF.


Kilograms (1000 g) - Select ON or OFF .


Avoirdupois Pounds (lb) - Select ON or OFF.

Ounces (oz) - Select ON or OFF.

Auto Zero Tracking - Scale will re-zero at 0.5d, 1d or 3d per second (0.5d is the LFT setting). Select 0.5d, 1d, 3d or OFF.



## 4.7 Print Menu

### Procedure

If you have just finished changing settings in the READ menu, press the **Data** button to go back to the READ MENU. READ should be blinking.

Press **M-** button to advance to the PRINT MENU. In this menu, you are able to change communication parameters and printing functions. Press the **Enter** button. Navigate through the menu as explained in section 4.3.1.


# Trooper Count Scale


## 4.7 Print Menu (Cont.)

### Procedure (Cont.)

Baud rate - Select either 1200, 2400, 4800, 9600 or 19200 baud.


Data Length - Select either 7 or 8.


Parity - Select either NONE, EVEN or ODD.


**NOTE:** No stop bit settings are required.


Print Part Number - Select ON or OFF


Print Count Quantity - Select ON or OFF.


Print APW - Select ON or OFF.


Print Net weight - Select ON or OFF.


Print Tare Value - Select ON or OFF.

Print Accumulated Count - Select ON or OFF.

Print Gross (Brutto) Weight - Select ON or OFF.


## 4.7 Print Menu (Cont.)

### Procedure (Cont.)

Select stability required for printing at 0.5d, 1d, 3d or OFF. OFF allows the printing when the weight is not stable. In LFT mode, stability is 1d.

Print Header - Allows printing of preset header that includes OHAUS, Model No. and unit serial number. Select ON or OFF.

Printout Style - See section 3.7.2 for examples. Select COL (24 column) or LINE (line).

## 4.8 Count Menu

### Procedure


If you have just finished changing settings in the Print menu, press the **Data** button to go back to the PRINT MENU. Press **M-** button to advance to the Count menu, then press **Enter** button. Navigate through the menu as explained in section 4.3.1.

Auto Add - Auto Add instructs the operator to add x pieces to reference count (see section 3.6). Select ON or OFF.

Print Clear - When Print clear is set to ON, the Scale prints the accumulated count memory and then clears the accumulated memory when the **PRINT Units** button is pressed. Select ON or OFF.

Auto Reference Optimization - If additional reference pieces less than or equal to the existing reference count are added to the pan, the Scale will recompute the APW based on the new reference weight (see section 3.5). Select ON or OFF.

Memory Change - A 10d weight change is required before the accumulated memory can be changed by pressing **M+** or **M-**. Select ON or OFF.


## 4.9 Lockout Menu


The Lockout menu (LOC) is a software-controlled feature that allows the Calibration, Setup, Read, Print and Count menus and the Library data to be locked to restrict tampering. When a menu or the Library is locked, the data can be viewed but not changed. To further restrict tampering, the Lock Switch on the circuit board can be set to prevent changes to the Lockout menu. See section 5.3 for Lock Switch information.

### Procedure


If you have just finished entering data in the Count menu, press the **Data** button to go back to the COUNT MENU. COUNT should be blinking. Press **M-** button to advance to the LOC menu. Press the **Enter** button. Navigate through the menu as explained in section 4.3.1.


LOC CAL - When set to ON, locks the CAL calibration sub-menu. Select ON or OFF.

LOC SETUP - When set to ON, locks the SETUP sub-menu. Select ON or OFF.


LOC READ - When set to ON, locks the READ sub-menu. Select ON or OFF.


LOC PRINT - When set to ON, locks the PRINT sub-menu. Select ON or OFF.

LOC COUNT - When set to ON, locks the COUNT sub-menu. Select ON or OFF.

LOC LIB - When set to ON, locks the Library records but still allows the accumulated memory to update (if activated). Select ON or OFF.

**NOTE:** Location 00 cannot be locked.

## 5. CALIBRATION AND SEALING

Although your Trooper Count scale has been calibrated at the factory, many conditions may require it to be re-calibrated when first put into use and periodically thereafter. Temperature change, geographic gravity variation, altitude change and abuse are a few reasons why a scale may need calibration.

*Span calibration* uses two points to adjust the scale. One point is called the zero value where there is no weight on the scale. The other point is the Span value that is often but not necessarily at full capacity. Trooper Count lets you choose between three Span values. These values are shown in the table below. The highest weight value generally gives the best performance over the entire weight range. Other span values may be chosen due to calibration weight availability or for improved accuracy at a lower weight range.

*Linearity calibration* uses three points to adjust the scale. In addition to the zero and full scale values an additional point is chosen between the two. The values are factory set and are shown in the table below. Linearity calibration should only be performed if determined to be necessary after a Span calibration.

*Before calibrating a scale* first determine that it is necessary. Place available calibration masses on the scale, as each mass is added the display should indicate the correct weight within +/- 1 scale division. If the scale is within this tolerance there is no need to calibrate the scale.

If the scale must be calibrated, use the table below to select the calibration value(s) that will be used. **Do not start the calibration procedure until you are certain that you have the appropriate calibration masses.**

Verify that the Cal menu is not locked (see section 4.9) and that the correct calibration unit (pound or kilogram) has been chosen in the Setup menu (See section 4.5).

**CALIBRATION MASSES**

<b>Cal in kg:</b>	<b>Span Cal choices</b>	<b>Linearity Cal (fixed)</b>
TC3RS	1, <b>2</b> , 3 kg	2 & 3 kg
TC6RS	2, <b>4</b> , 6 kg	4 & 6 kg
TC15RS	5, <b>10</b> , 15 kg	10 & 15 kg
TC30RS	10, <b>20</b> , 30 kg	20 & 30 kg
<b>Cal in lb:</b>	<b>Span Cal choices</b>	<b>Linearity Cal (fixed)</b>
TC3RS	2, <b>4</b> , 6 lb	4 & 6 lb
TC6RS	5, <b>10</b> , 15 lb	10 & 15 lb
TC15RS	10, <b>20</b> , 30 lb	20 & 30 lb
TC30RS	20, <b>40</b> , 60 lb	40 & 60 lb


Recommended tolerance Class for the calibration masses:  
ASTM Class 4 or OIML F2

# Trooper Count Scale


## 5.1 Span Calibration

### Procedure



The following example is for a 6kg scale.  
 With the Scale ON, press and hold the **SAMPLE Menu** button until CAL MENU is displayed.

Press the **Enter** button twice, NO should be blinking.

Press the **M-** button to change to YES, then press the **Enter** button.

Clear the pan, then press the **Enter** button.

(Scale is determining zero value)



Place the indicated calibration mass on the pan, then press the **Enter** button. The sample illustration indicates a 6kg Scale using a 2kg mass.

(Scale is determining span value)

If the calibration was successful, SPAN TEST DONE is displayed momentarily followed by CLEAR PAN THEN PRESS ENTER.

Clear the pan, then press the **Enter** button.


The scale will restart and return to weighing. If calibration was unsuccessful, refer to the troubleshooting section.


## 5.2 Linearity Calibration

### Procedure


The following example is with a 6kg Scale.  
With the Scale ON, press and hold the **SAMPLE Menu** button until CAL MENU is displayed.


Press the **Enter** button.



Press the **M-** button.


Press the **Enter** button.

Press the **M+** button, then press the **Enter** button.  
Follow display instructions.

(Scale is determining zero value)

Place the indicated calibration mass on the pan. Press the **Enter** button. The sample illustration indicates a 6kg Scale using a 4kg mass for the first linearity entry.

(Scale is determining first linearity value)


## 5.2 Linearity Calibration (Cont.)

### Procedure (Cont.)


Place the indicated calibration mass on the pan. Press the **Enter** button. The sample illustration indicates a 6kg Scale using a 6kg mass.


(Scale is determining second linearity value)

If the calibration was successful, CAL TEST DONE is displayed momentarily followed by CLEAR PAN THEN PRESS ENTER.

Clear the pan, then press the **Enter** button.

The scale will restart and return to weighing. If calibration was unsuccessful, refer to the troubleshooting section.

### 5.3 Legal for Trade Operation and Sealing

Legal for Trade (LFT) operation is possible through a hardware Lock Switch located on the PC board. The Scale setup parameters and calibration accuracy should be verified according to the local weights and measures standards prior to sealing.

#### Procedure

Set the LOC CAL, LOC Setup, LOC Read settings to ON per section 4.9.

Turn off the scale and disconnect the power. Remove the weighing platform.

Turn the Scale over in the position as shown and remove the Lock Switch Cover Plate.

Refer to the illustration at the right and notice the position of the LFT switch. To lock out the menus, slide the LFT switch to the position shown.

Replace the Lock Switch cover.

After the Scale has been tested and found to comply with applicable regulations, apply the appropriate seal.

#### LEAD AND WIRE SEAL


See illustration at right. Place wire seal through the holes in the security screws and ribs as shown and compress the lead seal in place.

#### PAPER SEAL


If paper seals will be used, place them over both access covers.

Turn the scale upright and replace the platform.


Reconnect power.

Bottom of Trooper Count

Sealing the Scale with a Lead Seal

Sealing the Scale with a Paper Seal

## 6. CARE AND MAINTENANCE

To keep the Scale operating properly, the housing should be kept clean and free from foreign material. If necessary, a cloth dampened with a mild detergent may be used.

### 6.1 Troubleshooting

If the Troubleshooting section does not resolve your problem, you need to contact an authorized Ohaus service agent.

SYMPTOM	PROBABLE CAUSE(S)	REMEDY
Unit will not turn on.	Adapter not plugged in or properly connected.	Check AC power.  Check AC Adapter connection.  Check connections at the back of the Scale.
	Batteries dead or not properly installed.	Check orientation of the batteries.  Replace batteries.
	Other failure.	Service required.
Cannot zero Scale, or will not zero when turned on.	Load on Scale exceeds allowable limits.	Remove load on Scale.
	Load cell damage.	Service required.
RS232 not working.	RS232 communication parameters not set correctly.	Verify communication parameters.
	Improper or loose cable connections.	Check cable connections.
Unable to calibrate unit.	Software Lockout switch set to ON.	Set Lock Cal to OFF in the Loc menu. Lock Switch on the circuit board may need to be set to off position.
	Incorrect value for calibration mass.	Use correct calibration mass.
Cannot display weight in desired weighing unit.	Desired unit not set to ON in Read menu.	Enable desired unit in Read menu.
Cannot change menu settings.	Menu has been locked.	Set selected menu to OFF in the LOC menu. Lock switch on the circuit board may need to be set to OFF position.
Key press causes three beeps and key press is ignored.	Improper key function or selection not available.	Review manual.

## 6.2 Error Codes List

The following list describes the various error codes and that can appear on the display.

- Low-BAT** [Batt ][Low ][ ] Indicates when batteries are weak.
- Error 1** [Err 1][over][Load] Indicates an overload condition.
- Error 2** [Err 2][Under][Load] Indicates an underload condition.
- Error 3** [Err 3][Low][APW] APW is below minimum allowed value.
- Error 4** [Err 4][Low][ref] Reference weight is below minimum value.
- Error 7** [Err 7][EEpr][Error] EEPROM data incorrect.
- Error 14** [Err 14][Zero][Error] Zero value exceeds the Zero % setting.
- Error 21** [Err 21][CAL][Error] Calibration incorrect.

## 6.3 Service Information

For Service assistance in the United States, please call Aftermarket, Ohaus Corporation, toll-free at (800) 526-0659. An Ohaus Product Service Specialist will be available to help you.

For Service assistance outside of the United states, call your local dealer.

## 6.4 Accessories

<u>Description</u>	<u>Ohaus Part No.</u>
RS232 Interface Cable/SF42 Printer	80500571
RS232 Interface Cable/PC 25 Pin	80500431
RS232 Interface Cable/PC 9 Pin	80500433
Printer	SF42

## 7. TECHNICAL DATA

### 7.1 Counting Limits

The following tables indicate the readability, minimum APW and minimum reference weight counting limits for the various models. Table 7-1 shows the limits with Auto-Add mode OFF. Table 7-2 shows the limits with the Auto-Add mode ON.

**TABLE 7-1. COUNTING LIMITS, AUTO ADD MODE OFF (1:100,000 Counting Resolution)**

Units = kg

Model	Readability	Min. APW	Min. Ref. Wt.
3kg	0.0005 kg	0.00003 kg	0.0005 kg
6 kg	0.001 kg	0.00006 kg	0.001 kg
15 kg	0.002 kg	0.00015 kg	0.002 kg
30 kg	0.005 kg	0.0003 kg	0.005 kg

Units g

Model	Readability	Min. APW	Min. Ref. Wt.
3 kg	0.5 g	0.03 g	0.5 g
6 kg	1 g	0.06 g	1 g
15 kg	2 g	0.15 g	2 g
30 kg	5 g	0.3 g	5 g

Units lb

Model	Readability	Min. APW	Min. Ref. Wt.
3 kg (6 lb)	0.001 lb	0.00006 lb	0.001 lb
6 kg (15 lb)	0.002 lb	0.00015 lb	0.002 lb
15 kg (30 lb)	0.005 lb	0.0003 lb	0.005 lb
30 kg (60 lb)	0.01 lb	0.0006 lb	0.01 lb

Units oz

Model	Readability	Min. APW	Min. Ref. Wt.
3 kg (96 oz)	0.02 oz	0.001 oz	0.02 oz
6 kg (240 oz)	0.05 oz	0.002 oz	0.05 oz
15 kg (480 oz)	0.1 oz	0.005 oz	0.1 oz
30 kg (960 oz)	0.2 oz	0.01 oz	0.2 oz

#### NOTES:

1. If **APW** is less than minimum, [Err3] [LOW] [APW] is displayed.
2. If **Ref. Wt.** is less than minimum, [Err4] [LOW] [REF] is displayed.

## 7.1 Counting Limits (Cont.)

TABLE 7-2. COUNTING LIMITS, AUTO ADD MODE ON (1:30,000 Counting Resolution)

Units = kg

Model	Readability	Min. APW	Min. Ref. Wt.	Auto-Add Min. wt.
3kg	0.0005 kg	0.0001 kg	0.001 kg	0.0025 kg
6 kg	0.001 kg	0.0002 kg	0.002 kg	0.005 kg
15 kg	0.002 kg	0.0005 kg	0.005 kg	0.010 kg
30 kg	0.005 kg	0.0010 kg	0.010 kg	0.025 kg

Units g

Model	Readability	Min. APW	Min. Ref. Wt.	Auto-Add Min. wt.
3 kg	0.5 g	0.1 g	1 g	2.5 g
6 kg	1 g	0.2 g	2 g	5 g
15 kg	2 g	0.5 g	5 g	10 g
30 kg	5 g	1 g	10 g	25 g

Units lb

Model	Readability	Min. APW	Min. Ref. Wt.	Auto-Add Min. wt.
3 kg (6 lb)	0.001 lb	0.0002 lb	0.002 lb	0.005 lb
6 kg (15 lb)	0.002 lb	0.0005 lb	0.005 lb	0.010 lb
15 kg (30 lb)	0.005 lb	0.001 lb	0.01 lb	0.025 lb
30 kg (60 lb)	0.01 lb	0.002 lb	0.02 lb	0.050 lb

Units oz

Model	Readability	Min. APW	Min. Ref. Wt.	Auto-Add Min. wt.
3 kg (96 oz)	0.02 oz	0.003 oz	0.03 oz	0.10 oz
6 kg (240 oz)	0.05 oz	0.010 oz	0.10 oz	0.25 oz
15 kg (480 oz)	0.1 oz	0.015 oz	0.15 oz	0.50 oz
30 kg (960 oz)	0.2 oz	0.030 oz	0.30 oz	0.10 oz

**NOTES:**

1. If **APW** is less than minimum, [Err3] [LOW] [APW] is displayed.
2. If **Ref. Wt.** is less than minimum, [Err4] [LOW] [REF] is displayed.
3. If **Auto-Add** weight is less than minimum, a calculated number of additional pieces are displayed [ADD] [xx] [PIECES].

# Trooper Count Scale

## 7.2 Specifications

Capacity x Readability	TC3RS	TC6RS	TC15RS	TC30RS
Enhanced	6 lb X 0.001 lb	15 lb X 0.002 lb	30 lb X 0.005 lb	60 lb X 0.01 lb
Enhanced	3 kg X 0.0005 kg	6 kg X 0.001 kg	15 kg X 0.002 kg	30 kg X 0.005 kg
Enhanced	3000 g X 0.5 g	6000 g X 1 g	15000 g X 2 g	30000 g X 5 g
Enhanced	96 oz X 0.02 oz	240 oz X 0.05 oz	480 oz X 0.1 oz	960 oz X 0.2 oz
Type Approved <sup>1</sup>	6 lb X 0.002 lb	15 lb X 0.005 lb	30 lb X 0.01 lb	60 lb X 0.02 lb
Type Approved <sup>1</sup>	3 kg X 0.001 kg	6 kg X 0.002 kg	15 kg X 0.005 kg	30 kg X 0.01 kg
Type Approved <sup>1</sup>	3000 g X 1 g	6000 g X 2 g	15000 g X 5 g	30000 g X 10 g
Type Approved <sup>1</sup>	96 oz X 0.05 oz	240 oz X 0.1 oz	480 oz X 0.2 oz	960 oz X 0.5 oz

### General Specifications

Platform Size (w x d) (in/mm)	13.5 X 9 / 340 X 230
Scale Dimensions (w x d x h) (in/mm)	13.5 X 14.5 X 4.5 / 343 X 368 X 114
Shipping Dimensions (w x d x h) (in/mm)	17 X 17 X 8.5 / 422 X 422 X 216
Weight (lb/kg)	8.6 / 3.9
Type Approved Resolution	1:3000
Maximum Resolution	1:6,000 - 1:7,500 depending on model
Stabilization time	1 - 4 seconds depending upon mode
Weighing units	lb, kg, g, oz
Display	3 backlit LCD Displays - 7 segment, 12 mm high
Power	AC Adapter or 6 alkaline C-type batteries
Typical Battery Life	Up to 100 hours
Span Calibration	1/3, 2/3, or full capacity
Linearity Calibration (3-point calibration)	0, 2/3 and full capacity
Auto-zero tracking	Off, 0.5, 1, or 3 divisions
Protection	IP42
Operating Temperature	41°F to 104°F / 5°C to 40°C
Storage Temperature	-40°F to 140°F / -40°C to 60°C
Construction:	
Internal Frame	Cast Aluminum
Weighing Platform	Stainless Steel
Housing	ABS
Keypad Overlay	Polyester

**Note:** All Trooper Count models are designed to meet the type approval requirements of NIST Handbook 44, Weights and Measures Regulations of Canada, OIML R76 and EEC EN45501.

### Ambient conditions

The technical data is valid under the following ambient conditions:

Ambient temperature	5°C to 40°C / 41°F to 104°F
Relative humidity	10%.....95%, noncondensing
Height above sea level	up to 4000m

## **LIMITED WARRANTY**

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.

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P/N 71141352 E © Ohaus Corporation 2003, all rights reserved / todos los derechos reservados / tous droits réservés.